

- ○ QC Wireless SDK Instructions for Use
 - ■ ■ update note:
 - 1.introduce
 - ■ 1.1The role of the SDK
 - 2. API description
 - 2.0 Access Conditions
 - 2.1 Permissions required by the SDK
 - 2.2 Access conditions
 - 2.3.1API
 - Scanning Devices
 - device connection:BleOperateManager.getInstance()
 - 2.3.2 Feature list:
 - Synchronize time, get the list of functions supported by the device
 - bracelet battery
 - Turn your wrist to brighten the screen
 - Do not disturb mode (command does not take effect if operating on the watch)
 - Time format, switch between public and English
 - weather
 - Continuous heart rate, blood oxygen, blood pressure switch
 - Set watch sports goals
 - Set the screen time of the watch, and customize the data type displayed by the watch face
 - Find equipment
 - Find phone features
 - Bracelet photo control
 - music control
 - Set the watch to factory reset
 - Sedentary reminder
 - Drink water reminder
 - Alarm data
 - firmware version number, hardware version number

- 2.3.3 Data synchronization:
 - Synchronize steps, distance, kcal for the day
 - Synchronized step data details
 - Sync Sleep Data Details
 - Synchronize the details of new sleep data and return according to SetTimeRsp
 - Sync heart rate data
 - Synchronized blood pressure function
 - Synchronized blood oxygen function
- 2.3.4 watch face
 - Read the list of watch faces in the watch
 - APP writes the dial to the watch
 - APP delete watch face
- 2.3.5 custom watch face
- 2.3.6 OTA upgrade function:
- 2.3.7 Manual measurement
- 2.3.8 One-click measurement

QC Wireless SDK Instructions for Use

1. Author: James
2. Shenzhen QC.wireless Technology Co., Ltd.
3. Version: 1.0.0

update note:

1. (2021/07/06) scan, connect, measure commands
2. (2021/07/20) Add setting command
3. (2021/07/21) Increase step count, heart rate, sleep data sync
4. (2022/02/28) Add new sleep algorithm

1.introduce

1.1The role of the SDK

Provide partner companies with the Android Bluetooth SDK for use with Green Orange wireless devices that provide basic and advanced functionality for a major watch or other device. This document is intended to explain the usage context, functionality, etc. of the API. Intended Audience and Reading Recommendations The intended audience and reader recommendations in this article are shown in Annex 1.

Reader	Role
Software Architecture Engineer	Architecture Analysis and Technical Guidance
Android development engineer	Have a certain android development ability, understand Ble related development technology

2. API description

2.0 Access Conditions

Android 5.0 or above, Bluetooth 4.0 or above.

2.1 Permissions required by the SDK

```

//network permissions
<uses-permission android:name="android.permission.INTERNET" />
//Bluetooth related permissions
<uses-permission android:name="android.permission.BLUETOOTH" />
<uses-permission
android:name="android.permission.BLUETOOTH_ADMIN" />
//Storage related permissions
<uses-permission
android:name="android.permission.WRITE_EXTERNAL_STORAGE" />
<uses-permission
android:name="android.permission.READ_EXTERNAL_STORAGE" />
//Location permission
<uses-permission
android:name="android.permission.ACCESS_COARSE_LOCATION" />
<uses-permission
android:name="android.permission.ACCESS_BACKGROUND_LOCATION" />
<uses-permission
android:name="android.permission.ACCESS_FINE_LOCATION" />

If it is android 12 or higher system
<uses-permission
android:name="android.permission.BLUETOOTH_CONNECT" />
<uses-permission
android:name="android.permission.BLUETOOTH_SCAN" />
<uses-permission
android:name="android.permission.BLUETOOTH_ADVERTISE" />

```

2.2 Access conditions

- Green Orange Wireless Wearables
- Green Orange Wireless SDK and documentation

2.3.1API

Scanning Devices

```

        //start scan
        BleScannerHelper.getInstance().scanDevice(final Context
context, UUID mUuid, final ScanWrapperCallback scanCallBack);
        //stop scan
        BleScannerHelper.getInstance().stopScan(Context context)
        //Specified device scan
        BleScannerHelper.getInstance().scanTheDevice(final Context
context, final String macAddress, final OnTheScanResult scanResult)

```

device connection:BleOperateManager.getInstance()

```

        //direct connection

        BleOperateManager.getInstance().connectDirectly(smartWatch.deviceAd
dress)
        //scan connections

        BleOperateManager.getInstance().connectWithScan(smartWatch.deviceAd
dress)
        //disconnect
        BleOperateManager.getInstance().unBindDevice()
        //reconnect
        BleOperateManager.getInstance().setNeedConnect(boolean
needConnect)
        //Called when bluetooth is turned off
        BleOperateManager.getInstance().setBluetoothTurnOff(false)
        BleOperateManager.getInstance().disconnect()
        //Turn on the system bluetooth monitor
        BleOperateManager.getInstance().setBluetoothTurnOff(true)

```

2.3.2 Feature list:

Synchronize time, get the list of functions supported by the device

```

        //set time
        CommandHandle.getInstance().executeReqCmd(SetTimeReq(0),
        ICommandResponse<SetTimeRsp>() {})
        //Callback Description
        public class SetTimeRsp extends BaseRspCmd {
            //support body temperature
            public boolean mSupportTemperature;
            //watch face
            public boolean mSupportPlate;
            //Support the menstrual cycle
            public boolean mSupportMenstruation;
            //Support custom watch faces
            public boolean mSupportCustomWallpaper;
            //Support blood oxygen
            public boolean mSupportBloodOxygen;
            //blood pressure support
            public boolean mSupportBloodPressure;
            //Support fatigue
            public boolean mSupportFeature;
            //Support one-key detection
            public boolean mSupportOneKeyCheck;
            //weather support
            public boolean mSupportWeather;
            //Support for new sleep protocol
            public boolean mNewSleepProtocol;
            //Supports up to 6 or 3 dials
            public int mMaxWatchFace;
        }

```

bracelet battery

```

        CommandHandle.getInstance().executeReqCmd(new
SimpleKeyReq(Constants.CMD_GET_DEVICE_ELECTRICITY_VALUE), new
ICommandResponse<BatteryRsp>() {

        @Override
        public void onDataResponse(BatteryRsp resultEntity) {
            if (resultEntity.getStatus() ==
BaseRspCmd.RESULT_OK) {
                //Get battery successfully
            }
        }

    });

    //Callback Description
    public class BatteryRsp extends BaseRspCmd {
        //battery [0-100]
        private int batteryValue;
    }

```

Turn your wrist to brighten the screen

```

    //Read the wrist roll setting in the watch
    CommandHandle.getInstance()
        .executeReqCmd(PalmScreenReq.getReadInstance(),
            ICommandResponse<PalmScreenRsp> {

                })

    //Write wrist roll settings isEnabled:true on, isLeft: true
    left hand
    CommandHandle.getInstance()
        .executeReqCmd(PalmScreenReq.getWriteInstance(boolean
isEnabled, boolean isLeft),
            ICommandResponse<PalmScreenRsp> {

                })

    PalmScreenReq Parameter Description
    /**
     * Whether the function is enabled
     */
    private boolean isEnabled;
    /**
     * Whether to wear on the left hand
     */
    private boolean isLeft;
    /**
     * whether to turn the wrist
     */
    private boolean needPalm;

```

Do not disturb mode (command does not take effect if operating on the watch)


```

//Read the watch into Do Not Disturb mode
    CommandHandle.getInstance()
        .executeReqCmd(DndReq.getReadInstance(),
            ICommandResponse<DndRsp> {

                })
    //Write Do Not Disturb

CommandHandle.getInstance().executeReqCmd(DndReq.getWriteInstance(b
oolean isEnabled,StartEndTimeEntity dndEntity,
    ICommandResponse<DndRsp> {

        })

StartEndTimeEntity Parameter Description
//start hour
private int startHour;
//start minute
private int startMinute;
//end hour
private int endHour;
//end minute
private int endMinute;

DndRsp Parameter Description
//function switch
private boolean isEnabled;
//Please refer to StartEndTimeEntity
private StartEndTimeEntity dndEntity;

```

Time format, switch between public and English

```

//Read watch parameters
CommandHandle.getInstance()
    .executeReqCmd(TimeFormatReq.getReadInstance(),
ICommandResponse<TimeFormatRsp> {})

//write to watch is24:true 24 false:12 metric:0 Metric
1: Imperial

CommandHandle.getInstance().executeReqCmd(TimeFormatReq.getWriteIns
tance(
    boolean is24, byte metric
),
    ICommandResponse<TimeFormatRsp> {

    })

TimeFormatRsp Parameter Description
//is24: true: 24 hours false:12hours
private boolean is24;
//isMetric: true: Metric false:Imperial
private boolean isMetric;

```

weather

```

//set weather units

CommandHandle.getInstance().executeReqCmd(DegreeSwitchReq.getWriteI
nstance(
    final boolean enable, final boolean isCelsius
),
    ICommandResponse<DegreeSwitchRsp> {

    }

)

//read weather units
CommandHandle.getInstance().executeReqCmd(DegreeSwitchReq
getReadInstance(),    ICommandResponse<DegreeSwitchRsp> {})

Parameter Description
enable: weather function true on, false off
isCelsius: true Celsius false Fahrenheit

```

```

    //Set the weather information, support the weather in the next 5
days

CommandHandle.getInstance().executeReqCmd(WeatherForecastReq.getWriteInstance(
    WeatherForecastBuilder builder),
ICommandResponse<WeatherForecastRsp>() { deviceResp ->
    if (deviceResp.status == BaseRspCmd.RESULT_OK) {

    }

    })

```

WeatherForecastBuilder Parameter Description

```

    //index Serial number 0-4 0 is today 1 is tomorrow ... support
up to 5 days
    private int index;
    //The unixtime of the weather is the second value at 0 o'clock
on a certain day
    private long timeStamp;
    //Weather type (0=unknown, 1=sunny, 2=cloudy, 3=rain, 4=snow,
5=haze, 6=lightning)
    private int weatherType;
    //Minimum temperature value (signed, range -128~127)
    private int minDegree;
    //Maximum temperature (with sign, range -128~127)
    private int maxDegree;
    //Humidity value (0~100)
    private int humidity;
    //Whether to bring an umbrella #01=with , others without
    private boolean takeUmbrella;

```

Continuous heart rate, blood oxygen, blood pressure switch

```

//Read continuous heart rate settings
CommandHandle.getInstance()
    .executeReqCmd(HeartRateSettingReq.getReadInstance(),
        ICommandResponse<HeartRateSettingRsp> {

        })
//Read Continuous SpO2 settings
CommandHandle.getInstance()
    .executeReqCmd(BloodOxygenSettingReq.getReadInstance(),
        ICommandResponse<BloodOxygenSettingRsp> {

        })
//Read continuous blood pressure settings
CommandHandle.getInstance()
    .executeReqCmd(BpSettingReq.getReadInstance(),
        ICommandResponse<BpSettingRsp> {

        })

//Write continuous heart rate switch isEnabled: true on,
false: off
CommandHandle.getInstance().executeReqCmd(
    HeartRateSettingReq.getWriteInstance(boolean isEnabled),
    ICommandResponse<HeartRateSettingRsp> {

    })

//Write continuous blood oxygen switch isEnabled: true on,
false: off
CommandHandle.getInstance().executeReqCmd(
    BloodOxygenSettingReq.getWriteInstance(boolean isEnabled),
    ICommandResponse<BloodOxygenSettingRsp> {

    })

//write blood pressure switch

CommandHandle.getInstance().executeReqCmd(BpSettingReq.getWriteInst
ance(boolean isEnabled, StartEndTimeEntity startEndTimeEntity, int
multiple), ICommandResponse<BpSettingRsp> {

    })
BpSettingRsp,Parameter Description
isEnabled: true on false off
StartEndTimeEntity The parameter description is the same as
above
multiple default 60

```

Set watch sports goals

```
CommandHandle.getInstance().executeReqCmd(  
    TargetSettingReq.getWriteInstance(  
        final int step, final int calorie, final int  
distance, final int sportMinute, final int sleepMinute  
    ), null  
)
```

Parameter Description

step: step target

calorie: Calorie target target, unit card, write kcal to a*1000

distance: Distance to target, in meters

sportMinute: Exercise duration target, in minutes

sleepMinute: Sleep duration target, in minutes

Set the screen time of the watch, and customize the data type displayed by the watch face

```
CommandHandle.getInstance().executeReqCmd(DisplayTimeReq.getWriteIn  
stance(final int displayTime, final int displayType, final int  
alpha),  
    ICommandResponse<DisplayTimeRsp> {  
  
    })
```

Parameter Description

displayTime Bright screen time range **[5,8]**

displayType Type of data displayed on the main **interface 0 steps 1
calories 2 weather 3 heart rate**

**alpha Screen brightness, temporarily useless, reserved, write 0 by
default**

Find equipment

```
CommandHandle.getInstance().executeReqCmd(FindDeviceReq(), null)
```

Find phone features

```
BleOperateManager.getInstance().addNotifyListener(  
    Constants.CMD_FIND_THE_PHONE.toInt(),  
    ICommandResponse<FindPhoneRsp>  
)  
FindPhoneRsp Parameter Description  
openOrClose 1: Start search 2: End search
```

Bracelet photo control

```

//The bracelet enters the camera interface

CommandHandle.getInstance().executeReqCmd(CameraReq(CameraReq.ACTION_
N_INT0_CAMARA_UI), null)
//The wristband is controlled by the bright screen on the camera
interface. The APP will send the bright screen command to keep the
watch bright. It is recommended to send it every 2 seconds.
    CommandHandle.getInstance().executeReqCmd(

CameraReq(CameraReq.ACTION_KEEP_SCREEN_ON),
            null
        )
//Bracelet click to take a photo event monitoring

BleOperateManager.getInstance().addNotifyListener(Constants.CMD_TAK
ING_PICTURE,new ICommandResponse<CameraNotifyRsp>(){

    @Override
    public void onDataResponse(CameraNotifyRsp
resultEntity) {

        }
    });
    resultEntity.getAction()
    Parameter Description
//The watch exits the camera interface
    CameraNotifyRsp.ACTION_FINISH
//The watch clicked on the photo event
    CameraNotifyRsp.ACTION_TAKE_PHOTO

CommandHandle.getInstance().executeReqCmd(CameraReq(CameraReq.ACTION_
N_FINISH), null)

```

music control

```

        //Send music information to the watch playing: true play, false
        pause progress play progress (default 0)
        //volume (current * 100 / max) name music song name
        FileHandle.getInstance().executeMusicSend(boolean playing, int
        progress, int volume, String name)

BleOperateManager.getInstance().addNotifyListener(Constants.CMD_MUS
IC_COMMAND,new ICommandResponse<MusicCommandRsp>(){

        @Override
        public void onDataResponse(MusicCommandRsp resultEntity)
    {

        }

    });
resultEntity.getAction() Parameter Description
1: play or pause
2: Previous song
3: Next song
4: Volume up
5: Volume down

```

- watch hang up

```

BleOperateManager.getInstance().addNotifyListener(Constants.CMD_PHO
NE_NOTIFY,new ICommandResponse<PhoneNotifyRsp>(){

        @Override
        public void onDataResponse(PhoneNotifyRsp resultEntity)
    {

        isReject() :true Hang up the phone
    }

    });

```

Set the watch to factory reset

```

CommandHandle.getInstance().executeReqCmd(RestoreKeyReq(Constants.C
MD_RE_STORE),null)

```

- message push


```
//The watch message push switch is fully turned on, and the APP  
should be actively opened  
    CommandHandle.getInstance().executeReqCmd(  
        SetANCSReq(), null  
    )  
//Send message push to watch  
    MessPushUtil.pushMsg(type,message:String)
```

PushMsgUintReq parameter description
type:

- 0x00: Call reminder 0x01: SMS reminder 0x02: QQ reminder
- 0x03: WeChat reminder,
- 0x04: incoming call to answer or hang up 0x05: Facebook
- message reminder 0x06: WhatsApp message reminder
- 0x07: Twitter message reminder 0x08: Skype message reminder
- 0x09: Line message reminder 0x0a: LinkedIn
- 0x0b: Instagram 0x0c: TIM message 0x0d: Snapchat
- 0x0e: others other types of notifications

Sedentary reminder

```

//Write a sedentary reminder

CommandHandle.getInstance().executeReqCmd(SetSitLongReq(StartEndTimeEntity startEndTimeEntity, byte weekMask, int cycle),
ICommandResponse<SimpleStatusRsp> { })
//Read sedentary reminders
CommandHandle.getInstance()

.executeReqCmd(SimpleKeyReq(Constants.CMD_GET_SIT_LONG),
                ICommandResponse<ReadSitLongRsp> { longSit ->

                })

```

StartEndTimeEntity parameter description is the same as above
weekMask repetition period
Bit0 = 0 means Sunday off, Bit0 = 1 means Sunday open.
Bit1 = 0 means closed on Monday, Bit1 = 1 means open on Monday.
Bit2 = 0 means close on Tuesday, Bit2 = 1 means open on Tuesday.
Bit3 = 0 means close on Wednesday, Bit3 = 1 means open on Wednesday.
Bit4 = 0 means close on Thursday, Bit4 = 1 means open on Thursday.
Bit5 = 0 means close on Friday, Bit5 = 1 means open on Friday.
Bit6 = 0 means closed on Saturday, Bit6 = 1 means open on Saturday.
The cycle reminder interval can only be one of 30, 60, and 90

Drink water reminder

```

    //Read the drinking water reminder index is the serial number of
the drinking water reminder [0, 8]

CommandHandle.getInstance().executeReqCmd(ReadDrinkAlarmReq(index),
ICommandResponse<ReadAlarmRsp>{})
    //write water reminder
    CommandHandle.getInstance().executeReqCmd(
        SetDrinkAlarmReq(AlarmEntity alarmEntity)
        null)
    }
AlarmEntity parameter description
    //Sequence number 0-8
    private int alarmIndex;
    //switch
    private int enable;
    // Reminder hour 0-23
    private int hour;
    //Reminder minute 0-59
    private int minute;
    // Repeat cycle as above
    private byte weekMask;

```

Alarm data

Note: Connect the watch to initialize:

LargeDataHandler.getInstance().initEnable() once, refer to the demo

//read the alarm clock

LargeDataHandler.getInstance().readAlarmWithCallback(**new**
IReadAlarmCallback() {}

//write the alarm clock

LargeDataHandler.getInstance().writeAlarm(entity);

AlarmNewEntity parameter description

//Total number of alarms

private int total;

//collection of alarm clocks

private List<AlarmBean> data;

AlarmBean parameter description

*//The length of the alarm clock: 4+content to byte[], refer to
demo*

private int alarmLength;

*//bit0~bit6 represents the switch corresponding to Monday to
Sunday, bit7 represents whether the current alarm clock is on*

private int repeatAndEnable;

//The number of minutes of the day, such as 480 represents 8:00

private int min;

//Alarm title, up to 30 bytes

private String content;

firmware version number, hardware version number

```

        //hardware information

CommandHandle.getInstance().execReadCmd(CommandHandle.getInstance()
    .getReadHwRequest());
        //firmware information

CommandHandle.getInstance().execReadCmd(CommandHandle.getInstance()
    .getReadFmRequest());

    Receiving implements this QCBluetoothCallbackCloneReceiver
    refer to demo MyBluetoothReceiver
    Judging UUID in the callback onCharacteristicRead

    override fun onCharacteristicRead(uuid: String?, data:
    ByteArray?) {
        if (uuid != null && data != null) {
            val version = String(data, Charsets.UTF_8)
            when(uuid){
                Constants.CHAR_FIRMWARE_REVISION.toString() -> {
                    //Firmware version number version
                }
                Constants.CHAR_HW_REVISION.toString() -> {
                    //hardware version number version number
version
                }
            }
        }
    }
}

```

2.3.3 Data synchronization:

Synchronize steps, distance, kcal for the day

```
CommandHandle.getInstance().executeReqCmd(  
    SimpleKeyReq(Constants.CMD_GET_STEP_TODAY),  
    ICommandResponse<TodaySportDataRsp> {})
```

TodaySportDataRsp parameter description

```
// days ago  
private int daysAgo;  
// date: year  
private int year;  
// date: month  
private int month;  
// date: day  
private int day;  
// total steps  
private int totalSteps;  
// running steps/aerobic steps  
private int runningSteps;  
// calorie value  
private int calorie;  
// walking distance  
private int walkDistance;  
// Movement time, in seconds  
private int sportDuration;  
// sleep time in seconds  
private int sleepDuration;
```

Synchronized step data details

```

        dayOffset Today 1: Yesterday 2: The day before yesterday,
supports synchronization for up to 7 days
        CommandHandle.getInstance().executeReqCmd(
            ReadDetailSportDataReq(dayOffset, 0, 95),
            ICommandResponse<ReadDetailSportDataRsp> {

                })
BleStepDetails parameter description
    //year
    private int year;
    //moon
    private int month;
    //day
    private int day;
    //15 minutes a point, the total number of points in a day is
    96 points, [0, 95], used to calculate the details of the number of
    steps per hour
    private int timeIndex=0;
    // calorie unit card
    private int calorie=0;
    //Step count
    private int walkSteps=0;
    //distance in meters
    private int distance=0;
    // keep for now
    private int runSteps=0;

```

Sync Sleep Data Details

```

    //deviceAddress Device mac address
    //dayOffset 0: Today 1: Yesterday 2: The day before yesterday,
    supports synchronization for up to 7 days
    //ISleepCallback callback

SleepAnalyzerUtils.getInstance().syncSleepReturnSleepDisplay(device
Address,dayOffset, ISleepCallback {
    //callback: SleepDisplay
})

SleepDisplay parameter description
    // total sleep time
    private int totalSleepDuration;
    // total time of deep sleep
    private int deepDuration;
    // total time of light sleep
    private int shallowDuration;
    // Go to sleep timestamp in seconds
    private int sleepTime;
    // wake up timestamp in seconds
    private int wakeTime;
    // A set of sleep data
    private List<SleepDataBean> list;
    private String address;

SleepDataBean parameter description
sleepStart start timestamp of a sleep in seconds
sleepEnd The end timestamp of a sleep in seconds
type sleep type 2: light sleep 1: deep sleep 3: awake

```

Synchronize the details of new sleep data and return according to SetTimeRsp


```

    //offset 0 today 1 yesterday
    public void syncSleepList(int offset, final
ILargeDataSleepResponse response)
    SleepNewProtoResp    parameter description
    //sleep start time
    private int st;
    //sleep end time
    private int et;
    //sleep collection
    private List<DetailBean> list;
    DetailBean parameter description
    // duration of a sleep type
    private int d;
    //type of sleep 2: light sleep 3: deep sleep 5: awake
    private int t;

```

Sync heart rate data

```

nowTime current time zone * 3600 + unix second value of current
time
    Sync yesterday: nowTime-(24*3600)*1,
    Sync the day before yesterday: nowTime-(24*3600)*2
    Data can be synchronized for up to three days
    val time = (getTimeZone() * 3600).toInt()
    val nowTime = date.unixTimestamp + time

    CommandHandle.getInstance().executeReqCmd(
        ReadHeartRateReq(nowTime),
        ICommandResponse<ReadHeartRateRsp> {

        })

ReadHeartRateRsp parameter description
    //nothing yet
    private int size = 0;
    //nothing yet
    private int index = 0;
    // unix second value of heart rate data
    private int mUtcTime;
    //The heart rate data array is one point every 5 minutes, the
    data subscript *5 is equal to the number of minutes of the day
    private byte[] mHeartRateArray;
    private boolean endFlag = false;

```

Synchronized blood pressure function

```
//Synchronized automatic blood pressure, measured once an hour
CommandHandle.getInstance()

.executeReqCmd(SimpleKeyReq(Constants.CMD_BP_TIMING_MONITOR_DATA),
ICommandResponse<BpDataRsp> {})
    BpDataEntity parameter description
    //year
    private int year;
    //month
    private int month;
    //day
    private int day;
    private int timeDelay;
    private ArrayList<BpValue> bpValues;

    BpValue parameter description
    //The minute of the day, usually the whole hour
    int timeMinute;
    //measured heart rate value
    int value;

    Get the blood pressure value calculated from the measured value
    //The heart rate value returned by the hr callback, age is the
    age of the user
    val sbp= CalcBloodPressureByHeart.cal_sbp(hr, age) (systolic
    blood pressure)
    //sbp heart rate calculated value
    val dbp=CalcBloodPressureByHeart.cal_dbp(sbp) (diastolic
    pressure)

    //Confirm blood pressure synchronization, call after receiving
    the callback, the watch will delete the records that have been
    synchronized

CommandHandle.getInstance().executeReqCmd(BpReadConformReq(true),nu
ll)

    //synchronize manual blood pressure
    CommandHandle.getInstance()
        .executeReqCmd(ReadPressureReq(0),
ICommandResponse<ReadBlePressureRsp> {})

    ReadBlePressureRsp.getValueList() parameter description
    BpPressure parameter description
    //time seconds value
    public long time;
```

```

public long time;
//(Diastolic pressure)

public int dbp;
//(systolic blood pressure)
public int sbp;

```

Synchronized blood oxygen function

```

LargeDataHandler.getInstance().syncBloodOxygenWithCallback(new
IBloodOxygenCallback() {
    @Override
    public void readBloodOxygen(List<BloodOxygenEntity>
data) {

        }
    });

BloodOxygenEntity parameter description
//data date
private String dateStr;
//Data minimum value array, one data per hour, a total of
24
private List<Integer> minArray;
//Data maximum value array, one data per hour, a total of
24
private List<Integer> maxArray;
//Data value at 0:00 on a certain day
private long unix_time;

```

2.3.4 watch face

Note: The watch supports up to 6 dial market dials and one custom dial. Among the 6 dial market dials, there are ones that can be deleted and those that cannot be deleted. The ones that can be deleted can be used to replace other dials. The dials cannot be added repeatedly. sex, can only add watch faces that do not exist in the watch

```

    //callback description: the callback of all operations on the
    dial
    inner class Callback : SimpleCallback() {
        //List of dial names returned by the device
        override fun onUpdatePlate(array: MutableList<PlateEntity>)
    {
        PlateEntity see description below
    }
        //The progress of the APP releasing the disk to the device
0-100
        override fun onProgress(percent: Int) {

        }
        //APP downloads the disk successfully
        override fun onComplete() {

        }
        //APP delete dial successfully
        override fun onDeletePlate() {

        }
        //Download the watch face or delete the watch face
exception
        override fun onUpdatePlateError(code: Int) {

        }
    }

    PlateEntity Description
    //Whether the dial can be deleted
    public boolean mDelete = false;
    //name of the dial
    public String mPlateName = "";

```

Read the list of watch faces in the watch

```

//Registry disk operation monitoring
FileHandle.getInstance().registerCallback(callback)
//initialization
FileHandle.getInstance().initRegister()
//Read the watch face list
FileHandle.getInstance().startObtainPlate()

```

APP writes the dial to the watch

```
// set file class
FileHandle.getInstance().currFileType =
FileHandle.TypeMarketWatchFace
//Registry disk operation monitoring
FileHandle.getInstance().registerCallback(callback)
//initialization
FileHandle.getInstance().initRegister()
//Execute file check, binPath is the path of the dial file
val prepare: Boolean =
FileHandle.getInstance().executeFilePrepare(binPath)
if (prepare) {
    //Start to issue the watch to the watch, the callback will
    return the progress and result
    FileHandle.getInstance().executeFileInit(name, 0x36)
}
```

APP delete watch face

```
// set file class
FileHandle.getInstance().currFileType =
FileHandle.TypeMarketWatchFace
//Registry disk operation monitoring
FileHandle.getInstance().registerCallback(callback)
//initialization
FileHandle.getInstance().initRegister()
//Delete the dial name The name of the dial, the watch will
return when reading
FileHandle.getInstance().executeFileDelete(name)
```

2.3.5 custom watch face

```
//Read custom dial coordinates, font color
LargeDataHandler.getInstance().readCustomizeDial(ILargeDataCallback
<CustomizeDialResp> respILargeDataCallback)
//CustomizeDialResp description
//The distance between the time and the left border is unit
pixels
private int timeLeft;
//The distance between time and upper border is unit pixel
private int timeTop;
```

```

    //The battery and date are unit pixels away from the left
border
    private int batteryLeft;
    //The unit pixel distance from the battery and date to the
upper border
    private int batteryTop;
    // The data and date are unit pixels away from the left border
    private int dataLeft;
    // The data and date are unit pixels away from the upper border
    private int dataTop;
    //font color
    private int r;
    private int g;
    private int b;

    //Write custom dial coordinates, font color
    LargeDataHandler.getInstance().writeCustomWatch(entity)
    CustomWatchFaceEntity parameter description Element type: 1 for
time 2: for battery and date 3: for data
    type: element type
    x: distance from the left border
    y: distance from the upper boundary
    r: font color R
    g: font color G
    b: font color B

    //Customize the dial to replace the back image
    FileHandler.getInstance().customizeWatchFace(Bitmap bitmap,int
width,int height,SimpleCallback callback)
    //Parameter Description
    Bitmap: background image
    width: watch width
    height: watch height
    SimpleCallback callback
    //Deliver the picture progress
    override fun onProgress(percent: Int) {

        }
    // Finished uploading the image
    override fun onComplete() {

        }

```

2.3.6 OTA upgrade function:

```

//dfu upgrade instance
val fuHandle= DfuHandle.getInstance()
//initialize callback
dfuHandle.initCallback()
//DFU file verification, path firmware file path
if (dfuHandle.checkFile(path)) {
    dfuHandle.start(dfuOpResult)
}
//dfuOpResult callback description
private val dfuOpResult: DfuHandle.IOpResult = object :
DfuHandle.IOpResult {
    override fun onActionResult(type: Int, errCode: Int) {
        if (errCode == DfuHandle.RSP_OK) {
            when (type) {
                1 -> dfuHandle.init()
                2 -> dfuHandle.sendPacket()
                3 -> dfuHandle.check()
                4 -> {
                    //The upgrade is successful, wait for the
device to restart
                    dfuHandle.endAndRelease()
                }
            }
        } else {
            //Upgrade exception or failure
        }
    }

    override fun onProgress(percent: Int) {
        // file upgrade progress
    }
}

```

2.3.7Manual measurement

```

        //StartHeartRateRsp parameter description
        private byte type; type 1: heart rate 2: blood
pressure 3: blood oxygen
        private byte errCode; measurement error code 0:
normal 1: measurement failed 2: measurement failed
        private byte value; measurement value: heart rate
or blood oxygen
        private byte sbp; blood pressure sbp
        private byte dbp; blood pressure dbp

        // manual heart rate
        BleOperateManager.getInstance().manualModeHeart(new
ICommandResponse<StartHeartRateRsp>() {
            @Override
            public void onDataResponse(StartHeartRateRsp
resultEntity) {

                }

        });
        //manual blood pressure
        BleOperateManager.getInstance().manualModeBP(new
ICommandResponse<StartHeartRateRsp>() {
            @Override
            public void onDataResponse(StartHeartRateRsp
resultEntity) {

                }

        });
        //manual blood oxygen
        BleOperateManager.getInstance().manualModeSpO2(new
ICommandResponse<StartHeartRateRsp>() {
            @Override
            public void onDataResponse(StartHeartRateRsp
resultEntity) {

                }

        });

```

2.3.8 One-click measurement


```
//Start the measurement (it is best to measure for more than 20  
seconds and then call the end, it is recommended 30 seconds)
```

```
lastFatigueValue The value of the last fatigue
```

```
lastFatigueTime Last test time in milliseconds
```

```
    BleOperateManager.getInstance().startOneKey(int  
lastFatigueValue,int lastFatigueTime,  
ICommandResponse<StartCalcDataRsp> callback)
```

```
//end measurement
```

```
BleOperateManager.getInstance().endOneKey()
```

StartCalcDataRsp parameter description

```
private byte errCode; 1 and 2 represent measurement failure
```

```
private byte heart; heart rate
```

```
private byte sbp; systolic pressure
```

```
private byte dbp; diastolic blood pressure
```

```
private int fatigue;
```

```
private int spo2; blood oxygen
```

```
private int score; score
```